


Anti-Chmp1B Rabbit Polyclonal Antibody For Wb, Elisa - Q7Lbr1

품목 번호: CM0025625



소개

CHMP1B Rabbit Polyclonal Antibody for Western blot (WB) and ELISA applications. Detects human, mouse, and rat CHMP1B, a key ESCRT-III component in endosomal sorting and cytokinesis.

[자세히 알아보기](#)

매개변수	값
제품명	CHMP1B Rabbit pAb
비고/별칭	Vps46B; C10orf2; C18orf2; CHMP1.5; Vps46-2; C18-ORF2; hVps46-2; CHMP1B
종(Species)	Human
유전자 ID (인간)	57132
유전자 ID	57132
면역원(Immunogen)	Recombinant fusion protein containing a sequence corresponding to amino acids 1-199 of human CHMP1B (NP_065145.2).
공급원(Source)	Rabbit
카테고리	Polyclonal Antibodies
애플리케이션(Application)	WB, ELISA
교차 반응성(Cross Reactivity)	Human, Mouse, Rat
SWISS	Q7LBR1
단백질 분자량	22kDa
배송(Shipping)	아이스백