

Monomethyl-Histone H3-K36 Rabbit Mab (Q16695/P68431) For Wb, Db, Ip, Elisa, Chip

품목 번호: CM0000361



소개

MonoMethyl-Histone H3-K36 Rabbit mAb (SWISS: Q16695/P68431) for epigenetics research. Validated in WB, DB, IP, ELISA, ChIP; human, mouse, rat reactivity. 15kDa, ice-bag shipping.

[자세히 알아보기](#)

매개변수	값
제품명	MonoMethyl-Histone H3-K36 Rabbit mAb
비고/별칭	H3/A; H3C2; H3C3; H3C4; H3C6; H3C7; H3C8; H3FA; H3C10; H3C11; H3C12; HIST1H3A; MonoMethyl-Histone H3-K36
종	인간
GeneID (인간)	8290/8350
GeneID	8290/8350
면역원	합성 펩타이드 인간 히스톤 H3 (NP_003520.1)의 K36 주변 합성 단일메틸화 펩타이드.
원천	토끼
범주	메틸화 항체
응용 분야	WB, DB, IP, ELISA, ChIP
교차 반응성	인간, 마우스, 래트, 기타 (광범위 예측)
SWISS	Q16695/P68431
단백질 무게	15kDa
배송	아이스백