

## Acetyl-Histone H3 K9 (Q16695/P68431) Rabbit Polyclonal Antibody - Wb Ip If Ihc Chip Elisa

품목 번호: CM0000170



### 소개

Acetyl-Histone H3-K9 rabbit polyclonal antibody for detection in human, mouse, rat. Suitable for WB, IP, IF/ICC, IHC-P, ChIP, ChIP-seq, ELISA. 15 kDa, broad cross-reactivity.

### [자세히 알아보기](#)

| 매개변수                     | 값                                                           |
|--------------------------|-------------------------------------------------------------|
| 제품명                      | Acetyl-Histone H3-K9 Rabbit pAb                             |
| 비고/별칭                    | H3t; H3.4; H3/g; H3FT; H3C16; HIST3H3; Acetyl-Histone H3-K9 |
| 종(Species)               | 사람(Human)                                                   |
| GeneID (사람)              | H3                                                          |
| GeneID                   | 8290/8350                                                   |
| 면역원(Immunogen)           | 사람 히스톤 H3 (NP_003520.1)의 아미노산 1-100 내 서열에 상응하는 합성 펩타이드      |
| 공급원(Source)              | 토끼(Rabbit)                                                  |
| 카테고리                     | 아세틸화 항체(Acetylated Antibodies)                              |
| 애플리케이션(Application)      | WB, IP, IF/ICC, IHC-P, ChIP, ChIP-seq, ELISA                |
| 교차 반응성(Cross Reactivity) | 사람(Human), 마우스(Mouse), 쥐(Rat), 기타 (광범위한 반응 예상)              |
| SWISS                    | Q16695/P68431                                               |
| 단백질 분자량                  | 15 kDa                                                      |
| 배송(Shipping)             | 아이스백                                                        |