

Anti-Cd13/Anpep 모노클로날 항체 (유세포분석용) - P15144

품목 번호: CM0007994



소개

rabbit monoclonal antibody against human CD13/ANPEP (P15144), validated for flow cytometry. Targets a broad-specificity aminopeptidase involved in peptide digestion, angiogenesis, and viral entry. Ideal for immunophenotyping and functional studies.

[자세히 알아보기](#)

파라미터	값
제품명	Rabbit anti-Human CD13/ANPEP mAb
비고/별칭	APN; AP-M; AP-N; CD13; LAP1; P150; PEPN; hAPN; GP150
종(Species)	Human
GeneID (Human)	290
GeneID	290
면역원(Immunogen)	Recombinant fusion protein containing a sequence corresponding to amino acids 69-967 of human CD13/ANPEP (NP_001141.2).
공급원(Source)	Rabbit
카테고리	모노클로날 항체
어플리케이션(Application)	FC
교차 반응성(Cross Reactivity)	Human
SWISS	P15144
단백질 분자량	110kDa
배송(Shipping)	아이스 팩(Ice bag)