

Anti-Cytokeratin 7 (Krt7) Rabbit Polyclonal Antibody For Wb, If-P, Elisa - P08729

품목 번호: CM0024480



소개

Cytokeratin 7 (KRT7) rabbit polyclonal antibody validated for WB, IF-P, ELISA. Reacts with human, mouse, rat. Targets a cytoskeletal protein involved in interferon signaling and viral translational regulation.

[자세히 알아보기](#)

매개변수	값
제품명	Cytokeratin 7 (KRT7) Rabbit pAb
비고/별칭	K7; CK7; SCL; K2C7; Cytokeratin 7 (KRT7)
종	Human
GeneID (인간)	3855
GeneID	3855
면역원	Recombinant fusion protein containing a sequence corresponding to amino acids 230-469 of human Cytokeratin 7 (KRT7) (NP_005547.3)
공급원	Rabbit
카테고리	Polyclonal Antibodies
애플리케이션	WB, IF-P, ELISA
교차 반응성	Human, Mouse, Rat
SWISS	P08729
단백질 분자량	51kDa
배송	Ice bag