


Anti-Human Pla2G7/Paf-Ah/Lp-Pla2 Rabbit Polyclonal Antibody For Wb And Elisa - Q13093

품목 번호: CM0023866



소개

Human PLA2G7/PAF-AH/Lp-PLA2 rabbit polyclonal antibody validated for WB and ELISA. Reacts with human and mouse. Useful for studying inflammatory phospholipid metabolism and cardiovascular disease research.

[자세히 알아보기](#)

매개변수	값
제품명	Human PLA2G7/PAF-AH/Lp-PLA2 Rabbit pAb
비고/별칭	PLA2G7; LDL-PLA2; LP-PLA2; PAFAD; PAFAH; phospholipase A2 group VII; Human PLA2G7/PAF-AH/Lp-PLA2
종	Human
GeneID (human)	7941
면역원	Recombinant fusion protein containing a sequence corresponding to amino acids 22-441 of human PLA2G7/PAF-AH/Lp-PLA2 (NP_005075.3)
원천	Rabbit
카테고리	Polyclonal Antibodies
응용 분야	WB, ELISA
교차 반응성	Human, Mouse
SWISS	Q13093
단백질 분자량	50kDa
배송	아이스백