


Anti-Mmp12 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - P39900

품목 번호: CM0023831



소개

MMP12 rabbit polyclonal antibody validated for WB, IF/ICC, ELISA in human, mouse, rat. Targets macrophage metalloelastase, involved in tissue remodeling and elastolysis. Ideal for inflammation and ECM studies.

[자세히 알아보기](#)

매개변수	값
제품명	MMP12 Rabbit pAb
비고/별칭	ME; HME; MME; MMP-12; MMP12
종(Species)	Human
유전자 ID (인간)	4321
유전자 ID	4321
면역원(Immunogen)	Synthetic peptide (amino acids 200-300 of human MMP12, NP_002417.2)
공급원(Source)	Rabbit
카테고리	Polyclonal Antibodies
애플리케이션	WB, IF/ICC, ELISA
교차 반응성	Human, Mouse, Rat
SWISS	P39900
단백질 분자량	54kDa
배송(Shipping)	Ice bag