


Anti-Nat10 Rabbit Monoclonal Antibody [Kd Validated] For Wb, Ip, If/Icc, Ihc-P, Elisa - Q9H0A0

품목 번호: CM0004693



소개

KD-validated rabbit monoclonal antibody targeting human NAT10 (RNA cytidine acetyltransferase). Suitable for WB, IP, IF/ICC, IHC-P, ELISA. Cross-reacts with mouse and rat. Ideal for studying RNA ac4C modifications and ribosome biogenesis.

[자세히 알아보기](#)

매개변수	값
제품명	Anti-NAT10 Rabbit Monoclonal Antibody [KD Validated]
비고/별칭	ALP; Kre33; NET43
종	사람(Human)
GeneID (사람)	55226
GeneID	55226
면역원	사람 NAT10 (NP_078938.2)의 아미노산 676-1025에 해당하는 서열을 포함하는 재조합 융합 단백질.
공급원	토끼(Rabbit)
카테고리	단일클론 항체(Monoclonal Antibodies)
애플리케이션	WB, IP, IF/ICC, IHC-P, ELISA
교차 반응성	사람(Human), 마우스(Mouse), 랫(Rat)
SWISS	Q9H0A0
단백질 분자량	116 kDa/107 kDa
배송	아이스백