

## Anti-Phospho-Akt1 (T450) Rabbit Monoclonal Antibody For Wb, Elisa - P31749

품목 번호: CM0000681



### 소개

Phospho-AKT1-T450 Rabbit mAb for WB and ELISA; specifically detects AKT1 phosphorylated at Thr450 in human and rat samples. Ideal for studying PI3K/AKT signaling in cancer, metabolism, and cell survival.

### [자세히 알아보기](#)

| 파라미터           | 값                                                             |
|----------------|---------------------------------------------------------------|
| 제품명            | Phospho-AKT1-T450 Rabbit mAb                                  |
| 비고/별칭          | AKT; PKB; RAC; PRKBA; PKB-ALPHA; RAC-ALPHA; Phospho-AKT1-T450 |
| 종              | 인간                                                            |
| GeneID (Human) | 207                                                           |
| GeneID         | 207                                                           |
| 면역원            | 인간 Akt1 (P31749)의 T450 주변 합성 인산화 펩타이드.                        |
| 원천             | 토끼                                                            |
| 카테고리           | 인산화 항체                                                        |
| 응용 분야          | WB, ELISA                                                     |
| 교차 반응성         | 인간, 쥐                                                         |
| SWISS          | P31749                                                        |
| 단백질 무게         | 56kDa                                                         |
| 배송             | 아이스백                                                          |