


Anti-Phospho-Hsp27/Hspb1-S78 Rabbit Monoclonal Antibody For Wb And Elisa - P04792

품목 번호: CM0000934



소개

Phospho-HSP27/HSPB1-S78 Rabbit Monoclonal Antibody validated for WB and ELISA. Specifically detects human HSPB1 phosphorylated at serine 78. Ideal for chaperone, stress response, and cancer research.

[자세히 알아보기](#)

파라미터	값
제품명	Phospho-HSP27/HSPB1-S78 Rabbit mAb
비고/별칭	CMT2F; HMN2B; HSP27; HSP28; Hsp25; SRP27; HS.76067; HEL-S-102; Phospho-HSP27/HSPB1-S78
종	인간
GeneID (인간)	3315
GeneID	3315
면역원	인간 HSPB1 (P04792)의 S78 주변 합성 인산화 펩타이드
원천	토끼
카테고리	인산화 항체
응용 분야	WB, ELISA
교차 반응성	인간
SWISS	P04792
단백질 분자량	23kDa
배송	아이스백