


Dimethyl-Histone H3-K9 Rabbit Polyclonal Antibody (Q16695/P68431) For Epigenetic Research | Wb, Chip, If/Icc

품목 번호: CM0000390



소개

DiMethyl-Histone H3-K9 rabbit polyclonal antibody (Q16695/P68431) validated for WB, ChIP, IF/ICC, and more. Reacts with human, mouse, rat, and other species. Ideal for epigenetic and chromatin studies.

[자세히 알아보기](#)

매개변수	값
제품명	DiMethyl-Histone H3-K9 Rabbit pAb
비고/별칭	H3t; H3.4; H3/g; H3FT; H3C16; HIST3H3; DiMethyl-Histone H3-K9
종(Species)	사람(Human)
GeneID (사람)	H3
GeneID	8290/8350
면역원(Immunogen)	합성 펩타이드: 사람 히스톤 H3 (NP_003520.1)의 아미노산 1-100 내 서열에 상응하는 합성 펩타이드
공급원(Source)	토끼(Rabbit)
카테고리	메틸화 항체(Methylated Antibodies)
응용(Application)	WB, DB, IHC-P, IF/ICC, ELISA, ChIP, ChIP-seq
교차 반응성(Cross Reactivity)	사람(Human), 마우스(Mouse), 쥐(Rat), 기타 (광범위한 반응 예상)
SWISS	Q16695/P68431
단백질 분자량	15 kDa
배송(Shipping)	아이스백(Ice bag)