


## Monomethyl-Histone H3-K27 (H3K27Me1) Rabbit Polyclonal Antibody (Q16695/P68431) For Wb, Ihc, Chip & Elisa

품목 번호: CM0000406



### 소개

Rabbit polyclonal antibody against monomethylated histone H3 at lysine 27 (H3K27me1). Validated for western blot, immunofluorescence, immunohistochemistry, dot blot, ChIP, and ELISA. Predicted cross-reactivity with human, mouse, and rat.

### [자세히 알아보기](#)

| 매개변수                     | 값                                                                                                        |
|--------------------------|----------------------------------------------------------------------------------------------------------|
| 제품명                      | MonoMethyl-Histone H3-K27 Rabbit pAb                                                                     |
| 비고/별칭                    | H3/A; H3C2; H3C3; H3C4; H3C6; H3C7; H3C8; H3FA; H3C10; H3C11; H3C12; HIST1H3A; MonoMethyl-Histone H3-K27 |
| 종(Species)               | 사람(Human)                                                                                                |
| GeneID (사람)              | H3                                                                                                       |
| GeneID                   | 8290/8350                                                                                                |
| 면역원(Immunogen)           | Synthetic peptide: A synthetic monomethylated peptide around K27 of human histone H3 (NP_003520.1)       |
| 공급원(Source)              | 토끼(Rabbit)                                                                                               |
| 카테고리                     | Methylated Antibodies                                                                                    |
| 애플리케이션(Application)      | WB, IF/ICC, IHC-P, DB, ChIP, ELISA                                                                       |
| 교차 반응성(Cross Reactivity) | Human, Mouse, Rat, Other (Wide Range Predicted)                                                          |
| SWISS                    | Q16695/P68431                                                                                            |
| 단백질 분자량(Protein Weight)  | 15 kDa                                                                                                   |
| 배송(Shipping)             | 아이스백(Ice bag)                                                                                            |