

Pe Mouse Anti-Human Cd56 Monoclonal Antibody For Flow Cytometry - P13591

품목 번호: CM0011175



소개

PE-conjugated mouse monoclonal antibody against human CD56 (NCAM1), validated for flow cytometry (FC). A reliable tool for studying natural killer cells, neural adhesion, and viral entry mechanisms. - P13591

[자세히 알아보기](#)

매개변수	값
제품명	PE Mouse anti-Human CD56 mAb
비고/별칭	CD56; NCAM; MSK39
종(Species)	Human
유전자 ID	4684
면역원(Immunogen)	합성 펩타이드 인간 CD56에서 유래한 합성 펩타이드
공급원(Source)	Mouse
카테고리	단클론 항체
응용(Application)	FC
교차 반응성	Human
SWISS	P13591
단백질 분자량	40kDa/73kDa/80kDa/83kDa/93kDa/95kDa
배송(Shipping)	아이스 팩