

Pe Rabbit Anti-Human Cd263/Trail-R3 Monoclonal Antibody For Fc - O14798

품목 번호: CM0004557

소개

PE-conjugated rabbit monoclonal antibody targeting human CD263/TRAIL-R3 (decoy receptor 1). Designed for flow cytometry (FC) to detect TRAIL-R3 on cell surfaces. Facilitates apoptosis and cancer research. - O14798

[자세히 알아보기](#)

|                          |                                                                                                                            |
|--------------------------|----------------------------------------------------------------------------------------------------------------------------|
| 매개변수                     | 값                                                                                                                          |
| 제품명                      | PE Rabbit anti-Human CD263/TRAIL-R3 mAb                                                                                    |
| 비고/별칭                    | LIT; DCRI; TRID; CD263; TRAILR3; TRAIL-R3; DCRI-TNFR                                                                       |
| 종(Species)               | Human                                                                                                                      |
| GeneID (human)           | 8794                                                                                                                       |
| GeneID                   | 8794                                                                                                                       |
| 면역원(Immunogen)           | Recombinant fusion protein containing a sequence corresponding to amino acids 26-235 of human CD263/TRAIL-R3(NP_003832.3). |
| 공급원(Source)              | Rabbit                                                                                                                     |
| 카테고리                     | Monoclonal Antibodies                                                                                                      |
| 애플리케이션(Application)      | FC                                                                                                                         |
| 교차 반응성(Cross Reactivity) | Human                                                                                                                      |
| SWISS                    | O14798                                                                                                                     |
| 단백질 분자량(Protein Weight)  | 27kDa                                                                                                                      |
| 배송(Shipping)             | Ice bag                                                                                                                    |