

Pe Rabbit Anti-Human Cd80/B7-1 Monoclonal Antibody For Fc - P33681

품목 번호: CM0005506



소개

PE-conjugated rabbit monoclonal antibody targeting human CD80 (B7-1), validated for flow cytometry. CD80 is a key costimulatory molecule on antigen-presenting cells critical for T-cell activation and immune checkpoint research.

[자세히 알아보기](#)

매개변수	사양
제품명	PE Rabbit anti-Human CD80/B7-1 mAb
비고/별칭	B7, BB1, B7-1, B7.1, LAB7, CD28LG, CD28LG1
종(Species)	Human
GeneID (인간)	941
GeneID	941
면역원(Immunogen)	인간 CD80/B7-1 (NP_005182.1)의 아미노산 35-242에 해당하는 서열을 포함하는 재조합 융합 단백질
공급원(Source)	Rabbit
카테고리	Monoclonal Antibodies
응용(Application)	FC
교차 반응성(Cross Reactivity)	Human
SWISS	P33681
단백질 분자량	33kDa
배송(Shipping)	아이스 백(Ice bag)