

Pe/Cyanine7 Anti-Human Cd27 Monoclonal Antibody For Flow Cytometry
- P26842

품목 번호: CM0008504



소개

PE/Cyanine7-conjugated rabbit monoclonal antibody against human CD27 (TNFRSF7), a costimulatory immune checkpoint receptor. Validated for flow cytometry (FC) applications. Ideal for T-cell and immune profiling research.

[자세히 알아보기](#)

매개변수	값
제품명	PE/Cyanine7 Rabbit anti-Human CD27 mAb
비고/별칭	T14; S152; Tp55; TNFRSF7; S152. LPFS2
종(Species)	Human
GeneID (인간)	939
GeneID	939
면역원(Immunogen)	재조합 단백질 인간 CD27 (NP_001233.2)의 아미노산 21-192에 해당하는 서열을 포함하는 재조합 융합 단백질.
공급원(Source)	Rabbit
카테고리	단클론 항체
응용(Application)	FC
교차 반응성(Cross Reactivity)	Human
SWISS	P26842
단백질 분자량	29kDa
배송(Shipping)	아이스 백(Ice bag)