


Phospho-A-Synuclein-Y133 Rabbit Pab For Wb, Elisa - P37840

품목 번호: CM0000937



소개

Phospho- α -Synuclein-Y133 Rabbit polyclonal antibody validated for WB and ELISA. Reacts with human and mouse samples. Useful for investigating synaptic activity, dopamine neurotransmission, and neurodegenerative disorders like Parkinson's and Alzheimer's disease.

[자세히 알아보기](#)

파라미터	값
제품명	Phospho- α -Synuclein-Y133 Rabbit pAb
비고/별칭	PDI; NACP; PARK1; PARK4; Phospho- α -Synuclein-Y133
종	인간
GeneID (인간)	6622
GeneID	6622
면역원	인간 알파-시뉴클레인 Y133 주변 잔기에 해당하는 합성 포스포-특이적 펩타이드
원천	토끼
카테고리	인산화 항체
응용 분야	WB, ELISA
교차 반응성	인간, 마우스
SWISS	P37840
단백질 분자량	14kDa
배송	아이스팩