


Phospho-Hdac4-S246/Hdac5-S259/Hdac9-S220 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - P56524 / Q9Uql6 / Q9Ukv0

품목 번호: CM0000428



소개

Phospho-specific rabbit polyclonal antibody detecting HDAC4 (S246), HDAC5 (S259), and HDAC9 (S220) for WB, IHC-P, ELISA. Reacts with human, mouse, rat. Suitable for chromatin regulation and transcription studies.

[자세히 알아보기](#)

매개변수	값
제품명	Phospho-HDAC4-S246/HDAC5-S259/HDAC9-S220 Rabbit pAb
비고/별칭	HDAC4/HDAC5/HDAC9; Phospho-HDAC4-S246/HDAC5-S259/HDAC9-S220
종	인간
GeneID (인간)	9759/10014/9734
GeneID	9759/10014/9734
면역원	합성 펩타이드 인간 HDAC4의 S246 주변 잔기에 해당하는 인산화 특이적 펩타이드
원산지	토끼
카테고리	인산화 항체
응용 분야	WB, IHC-P, ELISA
교차 반응성	인간, 마우스, 래트
SWISS	P56524/Q9UQL6/Q9UKV0
단백질 무게	57-65kDa/97-122kDa
배송	아이스백