

Phospho-Myosin Light Chain 2-T18 Rabbit Mab - P24844

품목 번호: CM0001915

소개

Phospho-Myosin Light Chain 2 (T18) Rabbit Monoclonal Antibody for WB and ELISA. Specific for human phosphorylated MYL9/MLC2, critical for smooth muscle and nonmuscle contractility research.

[자세히 알아보기](#)

매개변수	값
제품명	Phospho-Myosin Light Chain 2-T18 Rabbit mAb
비고/별칭	LC20; MLC2; MRLC1; MYRL2; MLC-2C; MMIHS4; Phospho-Myosin Light Chain 2-T18
종	인간
GeneID	10398
면역원	인간 Myosin Light Chain 2 (NP_006088.2)의 T18 주변 합성 인산화 펩타이드
원천	토끼
범주	인산화 항체
응용 분야	WB, ELISA
교차 반응성	인간
SWISS	P24844
단백질 무게	20kDa
배송	아이스백