

Rabbit Anti-Human Cd16B Mab For Fc - O75015

품목 번호: CM0008051



소개

-conjugated rabbit monoclonal antibody against human CD16b (FCGR3B). Designed for flow cytometry, detects low-affinity IgG Fc receptor on neutrophils. Useful for immunology and immune complex research.

[자세히 알아보기](#)

매개변수	값
제품명	Rabbit anti-Human CD16b mAb
비고/별칭	CD16, CD16-L, CD16b, FCG3, FCGR3, FCGR3B, FCR-10, FCRIII, FCRIIIb
종(Species)	Human
GeneID (Human)	2215
GeneID	2215
면역원(Immunogen)	Recombinant fusion protein containing a sequence corresponding to amino acids 20-208 of human CD16b (NP_000561.3)
공급원(Source)	Rabbit
카테고리	Monoclonal Antibodies
애플리케이션(Application)	FC
교차 반응성(Cross Reactivity)	Human
SWISS	O75015
단백질 분자량(Protein Weight)	26KDa
배송(Shipping)	Ice bag