

Rabbit Anti-Human/Mouse Integrin Alpha 5 (Itga5/Cd49E) Monoclonal Antibody For Flow Cytometry - P08648

품목 번호: CM0004140



소개

-conjugated rabbit monoclonal antibody against human/mouse integrin alpha-5 (ITGA5, CD49e), validated for flow cytometry. Functions as fibronectin receptor involved in cell adhesion. Suitable for IVD research and diagnostic assay development.

[자세히 알아보기](#)

매개변수	값
제품명	Rabbit anti-Human/Mouse Integrin alpha 5 (ITGA5/CD49e) mAb
비고/별칭	FNRA; CD49e; VLA-5; VLA5A
종(Species)	Human
유전자 ID (인간)	3678
유전자 ID	3678
면역원(Immunogen)	Recombinant protein of human/mouse alpha 5 (ITGA5/CD49e)
공급원(Source)	Rabbit
카테고리	Monoclonal Antibodies
응용(Application)	FC
교차 반응성(Cross Reactivity)	Human, Mouse
SWISS	P08648
단백질 분자량	115kDa
배송(Shipping)	Ice bag