

Rabbit Anti-Mouse Cd117/C-Kit Monoclonal Antibody For Flow Cytometry - P05532

품목 번호: CM0005778



소개

-conjugated rabbit monoclonal antibody against mouse CD117/c-Kit (SWISS: P05532). Suitable for flow cytometry (FC), recognizing the 22/108/109 kDa receptor. Ideal for stem cell and hematopoiesis research.

자세히 알아보기

사양	세부 정보
제품명	Rabbit anti-Mouse CD117/c-Kit mAb
비고/별칭	W; Bs; Fdc; Ssm; SC01; SC05; SOW3; CD117; c-KIT; Tr-kit; Gsfsc01; Gsfsc05; Gsfscow3
종	Mouse
유전자 ID (인간)	3815
유전자 ID	16590
면역원	마우스 CD117/c-Kit (NP_001116205.1)의 아미노산 25-523에 상응하는 서열을 포함하는 재조합 융합 단백질
공급원	Rabbit
카테고리	Monoclonal Antibodies
애플리케이션	FC
교차 반응성	Mouse
SWISS	P05532
단백질 분자량	22kDa/108kDa/109kDa
배송	아이스백