

Rabbit Anti-Mouse Cd6 Monoclonal Antibody For Flow Cytometry - Q61003

품목 번호: CM0009984

소개

-conjugated rabbit monoclonal antibody against mouse CD6 for flow cytometry. Recombinant immunogen (aa 20-155). Useful for investigating T-cell costimulation, immunological synapse maturation, and bacterial pattern recognition. - Q61003

[자세히 알아보기](#)

매개변수	값
제품명	Rabbit anti-Mouse CD6 mAb
비고/별칭	Eukaryota; Metazoa; Chordata; Craniata; Vertebrata; Euteleostomi; Mammalia; Eutheria; Euarchontoglires; Glires; Rodentia; Myomorpha; Muroidea; Muridae; Murinae; Mus; Mus
종(Species)	Mouse
GeneID (인간)	923
GeneID	12511
면역원(Immunogen)	Recombinant protein 마우스 CD6(NP_033982.3)의 아미노산 20-155에 해당하는 서열을 포함하는 재조합 융합 단백질.
공급원(Source)	Rabbit
카테고리	Monoclonal Antibodies
응용(Application)	FC
교차 반응성(Cross Reactivity)	Mouse
SWISS	Q61003
단백질 분자량	64kDa/68kDa/72kDa
배송(Shipping)	Ice bag